

ANCILLARY DEVICES

PIN-KAPS



You no longer have to use tape or tubing as our **BATRIK Pin-Kaps** will safely cover the sharp ends of your K-Wire and Steinmann Pins.

THEY ARE:

- Color coded by wire size which eliminates the selection of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Packaging designed to reduce waste
- Non-sterile (100 units/bag) also available
- Latex Free



ITEM NO. STERILE	ITEM NO. NON STERILE	COLOR	DIMENSIONS INCHES / MM	PACKAGING
14-1028	914-1028	Yellow	0.028" / 0.71 mm	40 pouches of 2
14-1035	914-1035	Blue	0.035" / 0.89 mm	40 pouches of 2
14-1045	914-1045	White	0.045" / 1.14 mm	40 pouches of 2
14-1062	914-1062	Green	0.062" / 1.57 mm	40 pouches of 2
14-1078	914-1078	Pink	0.078" / 1.98 mm	40 pouches of 2
14-1093	914-1093	Grey	0.094" / 2.36 mm	40 pouches of 2
14-1109	914-1109	Tan	0.109" / 2.77 mm	40 pouches of 2
14-1125	914-1125	Purple	0.125" / 3.17 mm	40 pouches of 2
14-1028-1	NA	Yellow	0.028" / 0.71 mm	40 pouches of 1
14-1035-1	NA	Blue	0.035" / 0.89 mm	40 pouches of 1
14-1045-1	NA	White	0.045" / 1.14 mm	40 pouches of 1
14-1062-1	NA	Green	0.062" / 1.57 mm	40 pouches of 1
14-1141	914-1141	Brown	0.078" / 1.98 mm	40 pouches of 2
14-1156	914-1156	Dark Blue	0.094" / 2.36 mm	40 pouches of 2
14-1188	914-1188	Red	0.109" / 2.77 mm	40 pouches of 2

SILICONE PIN MATS



Cushion, protect, and organize delicate surgical instruments for sterilization. Perforations allow for proper exposure and ventilation during sterilization. The Silicone Pin Mats can also be cut into smaller sections to fit any size basket or box.

- Perforated to allow proper exposure and ventilation during sterilization
- The user can easily cut the mat to fit their specific needs
- Color: Blue
- Latex Free

ITEM NO.	DIMENSIONS - INCHES / CM	PACKAGING
13-1000	13" x 18" / 33 x 45 cm	1 Unit



For NEW product launches and promotions visit www.batrik.com