

COLOR CODING SYSTEM

S.T.I.C.S.™ - SIMPLE TRACKING IDENTIFICATION COLOR SYSTEM

Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1.0" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2.0" (50.8 mm) x 44
- Total tape per sheet 374" (9500mm)

THEY ARE:

- Quickly identifies instruments by department or procedure
- Easy to apply
- Long lasting
- Tape adheres securely
- A free protective cover is provided with every sheet
- Latex Free

SOLID COLORS

1/4" (6.35 mm) Height / Various widths

ITEM NO.	COLOR
T-850-2-01	White
T-850-2-02	Royal Blue
T-850-2-03	Kelly Green
T-850-2-04	Red
T-850-2-05	Yellow
T-850-2-06	Purple
T-850-2-07	Orange
T-850-2-08	Brown
T-850-2-09	Black
T-850-2-10	Gray
T-850-2-11	Lavender
T-850-2-12	Beige
T-850-2-13	Light Blue
T-850-2-14	Lime Green
T-850-2-15	Pink
T-850-2-17	Plum
T-850-2-18	Mint Green

SPECIALTY STRIPE

1/4" (6.35 mm) Height / Various widths

ITEM NO.	COLOR
T-850-7-04	Yellow/Red

WHITE STRIPE

1/4" (6.35 mm) Height / Various widths

ITEM NO.	COLOR
T-850-4-02	White/Royal Blue
T-850-4-03	White/Kelly Green
T-850-4-04	White/Red
T-850-4-05	White/Yellow
T-850-4-07	White/Orange
T-850-4-08	White/Brown
T-850-4-11	White/Lavender
T-850-4-13	White/Light Blue
T-850-4-14	White/Lime Green
T-850-4-15	White/Pink
T-850-4-17	White/Plum

TRI-COLOR STRIPE

3/8" (9.53 mm) Height / Various widths

ITEM NO.	COLOR
T-850-5-51	Red/White/Blue
T-850-5-54	Red/Black/Green

BLACK STRIPE

1/4" (6.35 mm) Height / Various widths

ITEM NO.	COLOR
T-850-3-01	Black/White
T-850-3-03	Black/Kelly Green
T-850-3-04	Black/Red
T-850-3-05	Black/Yellow
T-850-3-07	Black/Orange
T-850-3-10	Black/Gray
T-850-3-11	Black/Lavender
T-850-3-12	Black/Beige
T-850-3-13	Black/Light Blue
T-850-3-14	Black/Lime Green
T-850-3-15	Black/Pink
T-850-3-17	Black/Plum

FOUR COLOR

1/4" (6.35 mm) Height / Various widths

ITEM NO.	COLOR
T-850-2-31	Red/Yellow/Purple/Orange
T-850-2-32	Blue/White/Black/Green

Custom requests are welcomed. Subject to volume and availability.
FOR MORE INFORMATION CALL 1.800.661.5432

IDENTIFICATION TAPE

BATRIK Identification Tapes are a simple and inexpensive identification solution for marking surgical instruments and withstands repeated exposure to sterilization.

- Quickly identify instruments by department or procedure
- Are easy to apply
- Withstands repeated exposure to sterilization
- Has more tape per roll (300" / 762 cm)
- Are available in 9 colors and 2 widths (1/4" / 0.635 cm and 1/8" / 0.3175 cm)
- Packaged in a convenient clear box.

Sterilization reports are available. All samples tested validated the efficiency of sterilization, even underneath the tape.



ITEM NO. 1/8" WIDE (0.3175CM)	ITEM NO. 1/4" WIDE (0.635CM)	COLOR
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-06	T-250-06	Purple
T-125-07	T-250-07	Orange
T-125-08	T-250-08	Brown
T-125-09	T-250-09	Black

