

IN-BRUSH™

TWISTED WIRE BRUSHES

TWISTED WIRE BRUSHES

Designed for use with suction tubes, our twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent brush cleaning capabilities.

NYLON TWISTED WIRE BRUSHES

We offer a wide selection of brushes for most of your instrument cleaning applications. Ideal for large bore, uniquely shaped, and endoscopic instruments.



OUR EXTENSIVE, COST EFFECTIVE PRODUCT LINE ALSO INCLUDES:

- Surgical Products
- Bite Blocks
- Ancillary Devices
- Cleaning Brushes
- Channel Brushes
- ...and much more!



FRAZIER / BARON SUCTION TUBES



ITEM NO.	SIZE	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
20-10079	6 FR	10" / 250 mm	1" / 25.4 mm	0.079" / 2.0 mm
20-10090	7 FR	10" / 250 mm	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	10" / 250 mm	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	10" / 250 mm	1" / 25.4 mm	0.118" / 3.0 mm
20-10129	10 FR	10" / 250 mm	1.5" / 38.0 mm	0.129" / 3.3 mm
20-10146	11 FR	10" / 250 mm	1.5" / 38.0 mm	0.146" / 3.7 mm
20-10157	12 FR	10" / 250 mm	1.5" / 38.0 mm	0.157" / 4.0 mm

PACKAGING: Bag of 3

NYLON TWISTED WIRE BRUSHES - 12" / 304 MM



ITEM NO.	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
22-12110	12" / 305 mm	2.0" / 51 mm	0.110" / 2.5 mm
22-12157	12" / 305 mm	3.0" / 76 mm	0.157" / 4.0 mm
22-12197	12" / 305 mm	3.0" / 76 mm	0.197" / 5.0 mm
22-12236	12" / 305 mm	3.0" / 76 mm	0.236" / 6.0 mm
22-12393	12" / 305 mm	3.0" / 76 mm	0.393" / 10.0 mm

PACKAGING: Bag of 3

NYLON TWISTED WIRE BRUSHES - 16" / 406 MM



ITEM NO.	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
23-16197	16" / 406 mm	2.50" / 63.5 mm	0.197" / 5.0 mm
23-16236	16" / 406 mm	2.50" / 63.5 mm	0.236" / 6.0 mm
23-16276	16" / 406 mm	3.00" / 76.0 mm	0.276" / 7.0 mm
23-16394	16" / 406 mm	3.00" / 76.0 mm	0.394" / 10.0 mm
23-16472	16" / 406 mm	3.00" / 76.0 mm	0.472" / 12.0 mm
23-16591	16" / 406 mm	3.00" / 76.0 mm	0.591" / 15.0 mm
23-16787	16" / 406 mm	3.00" / 76.0 mm	0.787" / 20.0 mm

PACKAGING: Bag of 3

NYLON TWISTED WIRE BRUSHES - 24" / 610 MM



ITEM NO.	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
24-24098	24" / 610 mm	2.50" / 63.5 mm	0.098" / 2.5 mm
24-24197	24" / 610 mm	2.50" / 63.5 mm	0.197" / 5.0 mm
24-24276	24" / 610 mm	3.25" / 82.6 mm	0.276" / 7.0 mm
24-24394	24" / 610 mm	3.25" / 82.6 mm	0.394" / 10.0 mm
24-24472	24" / 610 mm	3.25" / 82.6 mm	0.472" / 12.0 mm

PACKAGING: Bag of 3

BRONCHO BRUSHES



ITEM NO.	DESCRIPTION	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3

TRACHEOSTOMY BRUSHES



ITEM NO.	DESCRIPTION	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
25-08500	Large 6 - 12.3 3/4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000 - 5.1 7/8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3