

IN-BRUSH™

ALL PURPOSE BRUSHES AND FACE SHIELDS

PLASTIC HANDLE BRUSHES

We offer a wide selection of brushes for most of your instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.



PLASTIC SINGLE END NYLON - HANDLE BRUSH

ITEM NO.	DESCRIPTION	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
28-72550	Single End Nylon	7.25" / 184 mm	1.50" / 38 mm	0.50" / 12.7 mm

PACKAGING: Bag of 3 - ITEM NO.: 28-72550-50 / Bag of 50



PLASTIC SINGLE END STAINLESS STEEL - HANDLE BRUSH

ITEM NO.	DESCRIPTION	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
29-72550	Single End Stainless Steel	7.25" / 184 mm	1.50" / 38 mm	0.50" / 12.7 mm

PACKAGING: Bag of 3 - ITEM NO.: 29-72550-50 / Bag of 50



PLASTIC DOUBLE END NYLON - HANDLE BRUSH

ITEM NO.	BRUSH	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
28-69394	Large Brush-End	6.90" / 175 mm	1.38" / 35 mm	0.39" / 10.0 mm
	Small Brush-End	6.90" / 175 mm	1.00" / 25 mm	0.20" / 5.0 mm

PACKAGING: Bag of 3 - ITEM NO.: 28-69394-50 / Bag of 50

LARGE BOTTLE JAR CLEANING BRUSH



ITEM NO.	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
50-10295	10.40" / 264 mm	4.92" / 125 mm	2.95" / 75.0 mm

PACKAGING: Bag of 3

FACE SHIELD

- Flat viewing area with side closures which provides all-around greater protection.
- High Quality 6 mil crystal clear anti-fog Mylar, minimizes distortion, glare and provides excellent visibility.
- Comfortable self-shaping foam head-piece molds to the forehead.
- Other styles available.



ITEM NO.	PACKAGING
FI-160-100C	100 units / Box

IN-BRUSH™ PLASTIC HANDLE

Choose either the regular toothbrush style (nylon or stainless) or try our double-ended brush. Order in volume and save even more.



FLAT BRUSHES



ITEM NO.	BRISTLE TYPE	BRUSH LENGTH INCHES / MM	BRUSH L X I INCHES / MM
50-40915	Semi Rigid	6.692" / 170 mm	2.362"x0.787" / 60x20 mm
50-40916	Rigid	6.692" / 170 mm	2.362"x0.787" / 60x20 mm

PACKAGING: Bag of 3

HAND SCRUB NYLON BRUSH AUTOCLAVABLE



ITEM NO.	LENGTH INCHES / MM	BRUSH LENGTH INCHES / MM	BRUSH (Ø) INCHES / MM
50-40911	4.09" / 104 mm	3.74" / 95 mm	1.10" / 28.0 mm

PACKAGING: Bag of 3