

INSTRUMENT PROTECTION

INSTRUMENT PROTECTORS

OSTEOTOME PROTECTORS

- Specifically designed to protect during sterilization and handling
- Tips are protected from dulling prematurely or chipping
- Available in tinted blue perforated (holes)



ITEM NO.	DIMENSIONS INCHES	DIMENSIONS MM
OP-11	1/4" x 3/4"	6.35 mm x 19.05 mm
OP-12	3/8" x 1"	9.53 mm x 25.40 mm
OP-13	1/2" x 1"	12.70 mm x 25.40 mm
OP-14	3/4" x 1"	19.05 mm x 25.40 mm
OP-15	1" x 1"	25.40 mm x 25.40 mm
OP-16	1 1/4" x 1 1/4"	31.75 mm x 31.75 mm
OP-17	1 1/2" x 1 1/2"	38.10 mm x 38.10 mm
OP-18	1 3/4" x 1 1/4"	44.45 mm x 31.75 mm
OP-20	Assorted	

PACKAGING: 50 units per package - Non sterile



For NEW product launches and promotions visit www.batrik.com



RETRACTOR PROTECTORS SINGLE-TIP PERFORATED CLEAR "TUBE" STYLE

- Protect expensive retractors
- Steam and ETO compatible
- Latex Free

ITEM NO.	DIMENSIONS INCHES	DIMENSIONS CM	PACKAGING
RGS-01	1.5"	3.8 cm	55 units / bag
RGS-02	2.0"	5.1 cm	45 units / bag
RGS-03	3.0"	7.62 cm	40 units / bag

Non sterile

We pride ourselves on being extremely innovative, cost effective, and constantly updating our offerings to reflect the changing needs of our customers.

NEED TO DISCUSS A PRODUCT?
Call us today!
1.800.661.5432
(TOLL-FREE)

We're here to help you in any way that we can.

DOUBLE PROTECTORS

- Applied with the instruments in the open position thereby allowing the handles and jaws to be adequately exposed during the sterilization process
- Available in 6 tinted/perforated (holes)



ITEM NO.	DIMENSIONS INCHES	DIMENSIONS MM	COLOR
DG-01	1/16" x 3/4"	1.59 mm x 19.05 mm	White
DG-02	5/64" x 3/4"	1.98 mm x 19.05 mm	Blue
DG-03	7/64" x 3/4"	2.78 mm x 19.05 mm	Green
DG-04	1/8" x 1"	3.18 mm x 25.40 mm	Red
DG-05	3/16" x 1"	4.76 mm x 25.40 mm	Yellow
DG-07	1/16" x 3/8" x 1"	1.59 mm x 9.53 mm x 25.40 mm	Orange
DG-10	Assorted		

PACKAGING: 50 units per package - Non sterile
(Same value as purchasing 100 units)



RETRACTOR PROTECTORS "WEITLANER" PERFORATED "SQUARE" STYLE

- Protect expensive retractors
- Steam and ETO compatible
- Latex Free

ITEM NO.	SIZE	COLOR	PACKAGING
RGM-01	Small	Pink	50 units / bag
RGM-02	Small	Clear	50 units / bag
RGM-03	Medium	Green	50 units / bag
RGM-04	Medium	Clear	50 units / bag

Non sterile