




ATTENTION: PANDEMIC PRODUCT PREPAREDNESS COMMITTEE.
SUPERIOR FACE PROTECTION

Lock in low protected pricing
with guaranteed deferred delivery.



SF-75
SF-100



SF-180
SUPERIOR PROTECTION
Includes 6" Fluid Drape Barrier




SAVEFACE™ Shields use patented technology and superior faceshield material (7 mil thickness Mylar) to ensure excellent visibility and comfort.

USE FEWER SHIELDS
Due to higher quality materials and improved comfort.

ALL SAVEFACE™ SHIELDS HAVE:

- Patented flat viewing area that reduces distortion and glare.
- Self-shaping foam head piece comfortably molds to forehead.
- High quality 7 mil thickness mylar shield will result in a longer period of use.
(Always dispose of shields that become contaminated)
- Treated with superior Anti-Fog.
- Wide elastic band offers added comfort.
- Meets OSHA Blood Borne Pathogen Regulation.
- Latex Free

Product Ordering Information



100% Manufactured in North America

Serving you since 1993

Batrik Medical Manufacturing Inc.

Toll Free: 1.800.661.5432

Tel. no. 514.631.7988 • Fax: 514.631.9083 • Email: info@surgmed.com