

STERILISATION PRODUCTS

REUSABLE SYMMETRIC DOUBLE ARC™ WIRE BRUSHES (NOT ACTUAL SIZES)

ITEM No. 40-35012 - SAME SIZE BRUSH ON EACH END



12.5 MM

12.5 MM

14" / 35 CM SYMMETRIC DOUBLE WIRE BRUSHES

ITEM NO.	COLOR	LENGTH INCHES / CM	BRUSH (Ø) INCHES / MM	BRUSH LENGTH INCHES / MM	SHEATH (Ø) INCHES / MM
40-35005		12.59" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006		12.59" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008		12.59" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010		12.59" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012		12.59" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015		12.59" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020		12.59" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE WIRE BRUSHES

ITEM NO.	COLOR	LENGTH INCHES / CM	BRUSH (Ø) INCHES / MM	BRUSH LENGTH INCHES / MM	SHEATH (Ø) INCHES / MM
40-50005		19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006		19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008		19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010		19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012		19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015		19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020		19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE WIRE BRUSHES

ITEM NO.	COLOR	LENGTH INCHES / CM	BRUSH (Ø) INCHES / MM	BRUSH LENGTH INCHES / MM	SHEATH (Ø) INCHES / MM
40-65005		25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006		25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008		25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010		25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012		25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015		25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020		25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE WIRE BRUSHES

ITEM NO.	COLOR	LENGTH INCHES / CM	BRUSH (Ø) INCHES / MM	BRUSH LENGTH INCHES / MM	SHEATH (Ø) INCHES / MM
40-80005		31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006		31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008		31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010		31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80012		31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015		31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020		31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3