

SURGICAL PRODUCTS



GOLFF™ ANTI FOG SOLUTION

Laparoscopic surgical procedures often have problems with fogging. The **BATRIK Sterile Anti-fog Solution, GOLFF™**, is the solution that sets the standard for all other products.



GOLFF™:

- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Prevents the occurrence of fog on lenses
- Adhesive sponge for secure placement which is radio opaque
- Easy squeeze 5ml bottle
- OEM (available)
- Latex Free

ITEM NO.	DESCRIPTION	PACKAGING
GOLFF	Adhesive Sponge and Anti-fog per set	20 sterile sets / Box

CAUTERY TIP CLEANERS

Assists in removing material from cautery tips for both mono-polar and bi-polar.

THEY ARE:

- Foam and pad with abrasive material
- Adheres to disposable and non-disposable drapes
- Radio opaque
- Latex Free

DIMENSIONS
5 CM X 5 CM
2" X 2"



ITEM NO.	DESCRIPTION	PACKAGING
60-10000	Non Sterile - Cautery Tip Cleaners	2200 / Box
60-10100	Non Sterile - Cautery Tip Cleaners	1000 / Box
60-10200	Sterile - Cautery Tip Cleaners - Available in 2010	50 / Box



DID YOU KNOW...

BATRIK offers one of the largest selections of **ENDOSCOPY CHANNEL CLEANING BRUSHES** for surgical instruments! For more details... See page 12

SKINS™ - SURGICAL SKIN MARKERS

Designed with a beveled tip to enable the creation of multiple tracing lines (fine or bold) just by holding the pen at different angles.



THEY ARE:

- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Surgical prep resistant
- Latex Free

ITEM NO.	DESCRIPTION	PACKAGING
S-3000	STERILE - Multi-Tip Fine Marker with ruler and label set	50/box
S-3001	STERILE - Multi-Tip Fine Marker	50/box
S-93000	NON STERILE - Multi-Tip Fine Marker with ruler and label set	1000/box
S-93001	NON STERILE - Multi-Tip Fine Marker	1000/box
S-94000	NON STERILE - Blank ruler	100/bag
S-95000	NON STERILE - Set of 3 labels	100/bag
S-30005	STERILE - Multi-Tip Fine Marker with ruler and label set (Volume Packaging)	500/box
S-30015	STERILE - Multi-Tip Fine Marker (Volume Packaging)	500/box



SAVE by ordering our **SKINS™ MARKERS VOLUME PACKAGING** - more cost effective and reduces waste.

