

# ENDOSCOPY PRODUCTS

## REUSABLE EXTERNAL CLEANING BRUSH

### EXTERNAL CLEANING BRUSH FOR ENDOSCOPES/OPTICAL INSTRUMENTS

With its unique shape and innovative design, the **Batrik Reusable External Cleaning Brush** is the ultimate tool for cleaning the exterior of your endoscopic instruments without the risk of damaging expensive equipment:

- Easy-grip handle ensures comfort and control
- Numerous dense bristles provide additional surface contact
- Superior endoscopic scope cleaning with a brush specifically designed for this purpose

ITEM NO.	DESCRIPTION	PACKAGING
36-20010	Reusable External Cleaning Brush	3/Bag

HANDLE LENGTH DIAMETER  
120 MM/12 MM

BRUSH LENGTH/DIAMETER  
170 MM/ 25 MM

TOTAL LENGTH  
210 MM



- Polyamide (TYNEX®) bristles
- Superior bristle durability
- Steamable thermoplastic handle
- Ergonomic handle
- Steam/chemical sterilizable
- Latex Free

## SINGLE USE EXTERNAL CLEANING BRUSH

### EXTERNAL CLEANING BRUSH FOR ENDOSCOPES/OPTICAL INSTRUMENTS

An ideal fit for most endoscopic instruments, the **Batrik Single Use External Cleaning Brush** will provide excellent cleaning capabilities. Its "U" shaped bristle design provides continuous and effective cleaning. There are no sharp or metallic ends that can scratch or cause damage to valuable equipment.

Order in volume and **SAVE!**

ITEM NO.	DESCRIPTION	PACKAGING
36-20005	Single Use External Cleaning Brush	50/Box

HANDLE LENGTH DIAMETER  
80 MM/3 MM

BRUSH LENGTH/DIAMETER  
70 MM/ 25 MM

TOTAL LENGTH  
150 MM



- Many of the same features as the reusable brush
- Polypropylene bristles
- Ergonomic handle
- Vinyl caps
- Latex Free

### DID YOU KNOW...

Cleaning the exterior of your endoscope is an important step in the overall reprocessing of your endoscopic equipment.

### NEED ADDITIONAL INFORMATION?

Call us today!

**1.800.661.5432**

(TOLL-FREE)



Toll Free: 1.800.661.5432 Fax: 514.631.9083



850 Halpern Ave., Dorval (Montreal), Qc., Canada, H9P 1G6 • [www.batrik.com](http://www.batrik.com) • e-mail: [surgical@surgmed.com](mailto:surgical@surgmed.com)

# INSTRUMENT CLEANING BRUSHES

## REUSABLE PRECISION BRUSH

**EXTERNAL CLEANING BRUSH DESIGNED FOR CLEANING INSTRUMENT JAWS**

An ideal brush for the quick removal of difficult debris, the **Batrik Reusable Precision Brush** can also easily clean deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

**Eliminates the need for the stainless steel toothbrush-style brush, which can damage valuable instruments.**

ITEM NO.	DESCRIPTION	PACKAGING
50-40800	Reusable Precision External Cleaning Brush	3/Bag

HANDLE LENGTH DIAMETER  
120 MM / 12 MM



TOTAL LENGTH  
160 MM

BRUSH LENGTH / WIDTH / THICKNESS  
40 MM / 12 MM / 8 MM



- Polyamide (TYNEX®) bristles
- Bristles have an extra firm density and superior durability
- Steamable thermoplastic handle
- Ergonomic handle
- Steam/chemical sterilizable
- Latex Free



**CALL US TODAY TO RECEIVE YOUR FREE SAMPLE!**  
**1.800.661.5432 (TOLL-FREE)**

## REUSABLE JAR BRUSH

### REUSABLE JAR BRUSH

Practical and light, the **Batrik Reusable Jar Brush** has many long and flexible bristles to maximize large-surface cleaning. Its extra-long length also ensures a wider array of cleaning applications. The easy-grip handle also provides excellent manoeuvrability.

ITEM NO.	DESCRIPTION	PACKAGING
50-10290	Reusable Jar Brush	3/Bag

HANDLE LENGTH DIAMETER  
200 MM / 12 MM



TOTAL LENGTH  
330 MM

BRUSH LENGTH/DIAMETER  
110 MM / 75 MM



- Polyester derivative bristles
- Ergonomic handle
- Steamable
- Latex Free



Toll Free: 1.800.661.5432 Fax: 514.631.9083

850 Halpern Ave., Dorval (Montreal), Qc., Canada, H9P 1G6 • [www.batrik.com](http://www.batrik.com) • e-mail: [surgical@surgmed.com](mailto:surgical@surgmed.com)

