

COLOR CODING SYSTEM

PADLOCK SECURITY TAGS

An excellent product to secure instrument containers, security boxes, carts and cabinets from tampering. An additional advantage is that our printed padlocks allow you to identify instruments that require "sharpening" or "repair" thereby eliminating the re-processing of dull and broken instruments.

THEY ARE:

- Autoclavable* (*contact us for details*)
- Tamper-proof
- Available in 7 colors
- Seal must be broken to be removed



*Advise the temperature level so that we can recommend the appropriate tags.

NOTE: Tags can be consecutively numbered or custom printed



ITEM NO. SMALL SIZE	ITEM NO. MEDIUM SIZE	ITEM NO. LARGE SIZE	COLOR
15-501	15-601	15-551	White
15-502	15-602	15-552	Blue
15-503	15-603	15-553	Green
15-504	15-604	15-554	Red
15-505	15-605	15-555	Yellow
15-507	15-606	15-557	Orange
15-509	15-609	15-559	Black
	15-602ST		Blue "Sterile"
	15-604R		Red "Repair"
	15-605S		Yellow "Sharpen"
	15-605L		Yellow "Loaner"

PACKAGING: 100 per package



IMPORTANT INFORMATION:

This product range is being completely re-designed and re-engineered. Our new seals will withstand higher temperatures and accommodate a wide range of applications.

Target release date 02/2010.

For more information or details call
1.800.661.5432

INDICATOR SECURITY TAGS

Same features as our regular Padlock Security Tags, except these include an Indicator Dot which changes color to indicate exposure to either Steam or ETO sterilization.

FEATURES:

- Dual purpose indicator for ETO or Steam
- Indicator will change color based on the method of sterilization
- Available in 7 colors
- Large size will accommodate most instrument containers



NOTE: Tags can be consecutively numbered or custom printed



ITEM NO. MEDIUM SIZE	ITEM NO. LARGE SIZE	INCL. INDICATOR	COLOR
15-661	15-561	✓	White
15-662	15-562	✓	Blue
15-663	15-563	✓	Green
15-664	15-564	✓	Red
15-665	15-565	✓	Yellow
15-667	15-567	✓	Orange
15-669	15-569	✓	Black

PACKAGING: 100 per package



ACTUAL SECURITY TAGS MAY DIFFER FROM THE IMAGE. Contact us for details.